

## WUNDVERBÄNDE UND EIGENSCHAFTEN

| Produktmerkmale                      | Actolind® Hydrocolloid / Thin | Actolind® Hydrocolloid Foam | Actolind® Hydrocolloid Plus | Actolind® Plus Hydrocolloid | Actolind® Foam Non Adhesive | Actolind® Foam Adhesive | Actolind® Plus Foam | Actolind® Plus Foam Adhesive | Actolind® Wound Foam / Thin | Actolind® Wound Pad /Thin | Actolind® Wound Pad | Actolind® Wound Pad - Adhesive | Actolind® Wound Transparent Adhesive | Acilind® Sterile Gauze Compress | Acilind® Sterile Bandage |
|--------------------------------------|-------------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-------------------------|---------------------|------------------------------|-----------------------------|---------------------------|---------------------|--------------------------------|--------------------------------------|---------------------------------|--------------------------|
| CE 93/42/EWG Klassifizierung         | IIB                           | IIB                         | IIB                         | IIB                         | IIB                         | IIB                     | IIB                 | IIB                          | IIB                         | IIB                       | IIB                 | IS                             | IS                                   | IIA                             | IS                       |
| Einmalgebrauch, steril               | ✓                             | ✓                           | ✓                           | ✓                           | ✓                           | ✓                       | ✓                   | ✓                            | ✓                           | ✓                         | ✓                   | ✓                              | ✓                                    | ✓                               | ✓                        |
| Kleberand                            |                               |                             | ✓                           |                             |                             | ✓                       |                     |                              |                             |                           |                     | ✓                              | ✓                                    |                                 |                          |
| Leichte Anpassung an die Wunde       | ✓                             | ✓                           | ✓                           | ✓                           | ✓                           | ✓                       | ✓                   | ✓                            | ✓                           | ✓                         | ✓                   | ✓                              | ✓                                    | ✓                               | ✓                        |
| Saugfähigkeit                        | ✓ ❖                           | ✓ ❖                         | ✓ ❖                         | ✓ ❖                         | ✓ ❖                         | ✓ ❖                     | ✓ ❖                 | ✓ ❖                          | ✓ ❖                         | ✓ ❖                       |                     |                                |                                      | ✓ ❖ ❖                           |                          |
| Zum Zerschneiden geeignet            | ✓                             | ✓                           | ✓                           | ✓                           | ✓                           | ✓                       | ✓                   | ✓                            | ✓                           | ✓                         | ✓                   | ✓                              | ✓                                    | ✓                               | ✓                        |
| Durch Kontakt mit Exsudat Gelbildung | ✓                             | ✓                           | ✓                           | ✓                           | ✓                           | ✓                       |                     |                              | ✓                           | ✓                         |                     |                                |                                      |                                 |                          |
| Wasserdichte, atmungsaktive Schicht  | ✓                             | ✓                           | ✓                           | ✓                           | ✓                           | ✓                       | ✓                   | ✓                            | ✓                           | ✓                         |                     | ✓                              | ✓                                    |                                 |                          |
| Einzelteil-Ablösung                  | ✓                             | ✓                           | ✓                           | ✓                           | ✓                           | ✓                       | ✓                   | ✓                            |                             | ✓                         | ✓                   | ✓                              | ✓                                    | ✓                               | ✓                        |
| Selbstklebend                        | ✓                             | ✓                           | ✓                           | ✓                           |                             |                         |                     | ✓                            | ✓                           |                           |                     |                                |                                      |                                 |                          |
| Verschiedene Größen                  | ✓                             | ✓                           | ✓                           | ✓                           | ✓                           | ✓                       | ✓                   | ✓                            | ✓                           | ✓                         | ✓                   | ✓                              | ✓                                    | ✓                               | ✓                        |
| Dünne Option                         | ✓                             |                             |                             |                             |                             |                         |                     |                              |                             | ✓                         |                     |                                |                                      |                                 |                          |
| Barrierefunktion                     |                               |                             |                             |                             |                             |                         |                     |                              |                             |                           | ✓                   | ✓                              | ✓                                    |                                 | ✓                        |

❖ Hydrokolloid: Niedrig-Mittel  
❖ Foam: Multi-Medium  
❖ Gazekompress: Niedrig